# EM10-OWI-CHM-032 Rev. A

# **Ambient LOX Impact Testing**

NASA-STD-6001: Test 13A



Materials and Processes Laboratory Materials Test Branch, Building 4623

National Aeronautics and Space Administration George C. Marshall Space Flight Center Marshall Space Flight Center, AL 35812

Release Authority	Name	Title	Organization Date
Office of Primary Responsibility	[s] Gail H. Gordon	Materials Test Branch Chief	EM10 11/ <u>15/05</u>
	[s] Dennis Davis	Industrial Safety	QD50 11/ <u>16/05</u>



Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A Date: November 17, 2005 Page ii of iv			

Revision	Date	Originator	Description	Affected Pages
Baseline	2/4/05		Document converted from ED36-OWI-032. Previous history retained in system as part of canceled or superseded ISO Document files.	All
А	11/17/05		Hazardous Operations notification added to cover	Cover, ii

This document baselines the EM10 Organizational Work Instruction (OWI) for performing ambient-pressure liquid oxygen (LOX) impact tests in Marshall Space Flight Center's (MSFC's) Building 4623. Any deviation to this procedure shall be approved by the test engineer via an approved test plan. Any changes to the test equipment shall be noted on the tester maintenance log and approved by the test engineer. It is the responsibility of the test engineer to obtain NASA Contracting Officer's Technical Representative (COTR) approval where necessary for changes to the test equipment.

Any change to this OWI shall be submitted to and approved by the Materials Test Branch Chief, EM10. Revisions may also be submitted to the concurring organizations listed below for review and concurrence by memo. The original OWI and all changes shall be maintained by EM10. Any change to materials used requires a change to mechanical drawings, in addition to Chemistry Team Leader approval. All documentation shall be approved by the appropriate persons mentioned above and incorporated into the OWI before operation of the reconfigured test equipment can resume.

Concurring organizations:
Building 4623 Test Operations Contractor
Chemistry Team Leader
Environmental Health, AD60M

I	Marshall Space Flight Center Organizational Work Instruction			
ı	EM10			
ı	AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
ı	NASA-STD-6001: Test 13A Date: November 17, 2005 Page jij of jy			

# **Contents**

	Document	History Log	iii
1.0	1.1 So 1.2 Pu	ropeplicability	1 1
2.0	Applicable	Documents	2
3.0	3.1 D	efinitions cronyms	3
4.0	4.1 Si 4.2 P 4.3 Ei 4.4 Si 4.5 D 4.6 Si	sample Preparation	5 6 7 12
5.0	Notes		15
6.0	6.1 H 6.2 Si 6.3 Si ai 6.4 Ei 6.5 Ai 6.6 Ei	cautions and Warning Notes	16 17 18 19 19
7.0	7.1 A	ts, Data, Reports, and Formstachment: Striker Pin Boxing Guide	20
8.0	8.1 M 8.2 C	lemoranda alibration Records ecords	25 25
9.0	9.1 Si 9.2 Pi 9.3 R 9.4 Ci	pment, and Materialstandard Configuration of Tester	26 26 26
10.0	Personnel	Training	31
	Emergency	Phone Numbers	34

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page iv of iv	

This page is deliberately blank.

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A			

1.0 Scope

#### 1.1 Scope

The scope of this operating instruction is NASA-STD-6001, Test 13A, Mechanical Impact for Materials in Ambient Pressure LOX, as performed in Building 4634 at Marshall Space Flight Center.

## 1.2 Purpose

The purpose of Test 13A is to determine the mechanical impact sensitivity of materials when exposed to liquid oxygen at ambient pressure.

#### 1.3 Applicability

This instruction applies to the Chemistry Team, Materials Test Branch, of the Materials and Processes Laboratory.

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 2 of 34	

# 2.0 Applicable Documents

ASTM D 2512-82. Standard Test Method for Compatibility of Materials with Liquid Oxygen (Impact Sensitivity Threshold and Pass-Fail Techniques). 1982.

EM10-OWI-CHM-042. Test Sample Preparation for Testing in Building 4623.

EM10-OWI-CHM-050. Building 4623 Guidelines for Test Operations.

EM10-OWI-CHM-051. Receipt, Handling, Prioritizing, and Data Requirements of Samples Submitted for Testing in Building 4623 of the Materials and Processes Laboratory.

EM10-OWI-CHM-058. Chemical Hygiene Plan for Building 4623.

MPD 1840.3. MSFC Respiratory Protection Program.

MPR 1040.3. MSFC Emergency Plan.

MPR 1840.2. MSFC Hazard Communication Program.

MPR 8715.1. MSFC Safety, Health, and Environmental (SHE) Program.

MSFC-SPEC-164B. Specification for Cleanliness of Components for Use in Oxygen, Fuel, and Pneumatic Systems.

MWI 3410.1. Personnel Certification Program.

MWI 8621.1. Close Call and Mishap Reporting and Investigation Program.

NASA-STD-6001. Flammability, Odor, Offgassing, and Compatibility Requirements and Test Procedures for Materials in Environments That Support Combustion.

NHB-8060.1B. Flammability, Odor, Offgassing, and Compatibility Requirements and Test Procedures for Materials in Environments That Support Combustion.



**Note:** Always reference the current version of an applicable document.

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 3 of 34	

# 3.0 Definitions

#### 3.1 Definitions

Blank. A one-piece or two-piece specimen cup containing only an insert.

*Disk.* An anodized aluminum-alloy disk used as a carrier for dyes, dye penetrants, and emulsifiers.

Full access. A term meaning, "Anyone may enter the test cell."

*Insert.* A 300-series stainless-steel disk placed in the bottom of the specimen cup.

*Limited access*. A term meaning, "Only the test operator shall enter the test cell with appropriate personal protective equipment."

*Mechanical impact*. Energy delivered to a sample by a plummet that has been dropped from a pre-established height onto a striker pin in contact with the sample.

*NASA*. Marshall Space Flight Center EM10 responsible personnel.

No access. A term meaning, "No one shall enter the test cell."

One-piece insert cup (grease cup). An aluminum-alloy cup placed in the bottom of the one-piece specimen cup. Used for liquid and semisolid samples. Occasionally, stainless-steel insert cups are used.

Oxygen deficient. Oxygen level in air is less than 19.5%.

Oxygen enriched. Oxygen level in air is greater than 23.5%.

*Reaction.* A chemical change or transformation in the sample resulting from mechanical impact and determined by an audible report, a visually detected flash, sustained burning, or obvious charring of the sample, sample cup, or striker pin.

Specimen cup. An aluminum-alloy cup used for solid and liquid samples.

*Tag out.* Placement of a tag-out device on an energy-isolating device to indicate that the energy-isolating device and equipment being controlled shall not be operated until the tag-out device is removed by the person who placed it there.

*Test 13A.* Mechanical impact testing for materials in ambient pressure LOX per NASA-STD-6001.

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 4 of 34	

*Test area.* The portion of Building 4623 and fenced area south of the north wall of Room 122.

*Test cell.* Room 122A of Building 4623, which contains the impact testing equipment.

*Test engineer*. The person responsible for correctly following the approved test plan for a specific test -- from sample receipt to test data evaluation.

*Test operator*. The person responsible for conducting a test under the guidance of the test engineer.

#### 3.2 Acronyms

$GN_2$	Gaseous nitrogen
LED	Light Emitting Diode
$LN_2$	Liquid nitrogen
LOX	Liquid oxygen
MSDS	Material Safety Data Sheet
<i>MSFC</i>	Marshall Space Flight Center
OWI	Organizational Work Instruction
PTFE	Polytetrafluoroethylene

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A	
NASA-STD-6001: Test 13A			

## 4.0 Instructions

All operations of this equipment shall be conducted using the applicable documents referenced above (section 2). All data and test results shall be recorded on form EM10-F-CHM-013, the Test 13A Test Data Sheet (section 7.2, Figure 2). A summary of pertinent test information and test results shall be compiled in a memorandum, signed by the test organization management, and mailed to the test requester.

#### 4.1 Sample Preparation

The sample preparation technician shall prepare Test 13A samples according to EM10-OWI-CHM-042, Test Sample Preparation for Testing in Building 4623. When non-standard samples are received, the sample preparation technician shall follow the directions written in the test plan for that test request. If this information is not provided with the test plan, the sample preparation technician shall seek clarification from the test engineer.

Before testing begins, the *test operator* **shall review** the information supplied on the test data sheet (prepared by the sample preparation technician) to make certain the information is complete and appears sound. *If a problem is identified*, the *test operator* **shall notify** the test engineer. The *test operator* shall also:

- **Verify** that the test request number and material designation are identical on all paperwork.
- Confirm that the prepared samples agree with the test request.
- **Verify** that the sample preparation technician has noted if the sample has been cleaned or if the sample does not need to be cleaned.
- **Note** any flaws or imperfections in the sample, and **record** these on the test data sheet.
- **Review** the signed test plan and the original test request before proceeding. *If the test plan and the test request do not agree*, **request** clarification from the test engineer.

#### 4.2 Pre-Test Photography

The sample preparation technician **shall take** a pre-test photograph of at least one of the samples and **place** three copies of the photograph in the test folder. If the pre-test photograph has not been taken, the test operator **shall take** the photograph and **place** three copies of the photograph in the test folder before proceeding with the test. The entire sample shall be visible in the photo. Steps for photographing samples are outlined in the Photography Operating Guide.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 6 of 34

#### 4.3 Equipment Checkout

At the beginning of the test day, the *test operator* **shall perform** the following steps:

- 4.3.1. **Observe** at the oxygen monitor station located in Room 116. **Locate** the monitor for Room 122A. **Ensure** that the digital LED readout indicates a normal oxygen level (19.5 to 23.5%) before entering the test cell.
- 4.3.2. **Ensure** that the test cell floor and the room behind the test cell are visibly clean. *If conditions warrant*, **scrub** the floor of the cell with a detergent solution, and **rinse** with water.
- 4.3.3. Check the plummet assembly before each test request. Ensure that the plummet spiders are securely attached to the plummet body and that bearing wheels are securely mounted on the spiders; use a wrench to tighten these as needed. If any bearings are missing, replace these before testing. Ensure that the spiders are in proper alignment with the test stand columns and that bearing wheels rotate freely. Verify that the clamp on the magnet housing is secure. (Use a wrench if necessary.) Ensure that all guide track retaining bolts are hand-tight. If not, notify the test engineer who shall call the Fabrication Group to come align the track(s) and tighten the bolt(s). To ascertain the alignment of the vertical guide tracks and the error caused by friction of the plummet assembly, time the plummet fall from a given height. For this test, the free-fall timing is to be within 3% of the value given for the associated height (section 7.2, Table 1).
- <u>4.3.4.</u> **Check** the anvil plate for indentations or other irregularities. **Examine** the concrete pedestal for any crumbling. **Examine** the rails for bends or large scratches. **Report** any problems to the test engineer.
- <u>4.3.5.</u> **Ensure** the proximity switches (timing sensors), electromagnet, mechanical safety catch, exhaust fan, and the facility warning beacon are working. **Report** problems to the test engineer before making any adjustments.
  - **Refer** to section 7.2, Table 1, for drop times.
  - Raise the plummet to the magnet base to test the electromagnet and safety catch.
  - **Turn on** the exhaust fan, and **listen** for operation.
  - **Turn on** the facility warning beacon.



**CAUTION:** Verify proper operation of the safety catch before working in the plummet fall/impact path.

<u>4.3.6.</u> **Thoroughly clean** the guide tracks, plummet, anvil plate, striker-pin guide, specimen-cup holder, baseplate, sample holder top plate, hand dewars, forceps, sample freezing box, and measuring rods to be used. **Use** cheesecloth soaked with deionized water and non-ionic soap. **Rinse** with deionized water. *If the* 

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page <b>7</b> of <b>34</b>

equipment is oily or greasy, **clean** it again with fresh trichlorotrifluoroethane. Trichloroethylene or Brulin<sup>®</sup> 815 MX may be substituted for trichlorotrifluoroethane, *if the cleaning is performed in a well-ventilated area*. **Allow** parts to dry completely before assembly.





**Note**: **Discard** any rags and debris soaked with trichlorotrifluoroethane (or organic materials) in the *Organic Rags and Debris* hazardous control drum, which is located in Room 128 (in the fenced area south) of Building 4623. **Do not discard these rags in any other receptacle.** 



4.3.7. Clean the pre-cooling moat by placing all parts in the ultrasonic cleaner with a mixture of non-ionic soap and deionized water for 20 minutes. Rinse with deionized water. If parts are still dirty or need further cleaning, wipe each part with trichlorotrifluoroethane. Trichloroethylene or Brulin® 815 MX may be substituted for trichlorotrifluoroethane, provided the cleaning is performed in a well-ventilated area. Allow parts to dry completely before assembly. Clean and store equipment after use.





### 4.4 System Setup and Sample Loading


The test operator **shall perform** the following steps:

- 4.4.1. **Inspect** the contents of the current test folder to ensure that the following information is provided:
  - Test request
  - Signed test plan
  - Test data sheets
  - Sample preparation sheet
  - Test material's MSDS or the Exclusion Statement for the material/component being tested
  - Pre-test photographs (**See** section 4.2.)
  - Ambient-Pressure LOX Impact Tester Pre-Test Checklist.

WARNING: Read the test material's MSDS to ensure familiarity with all safety precautions associated with the material. Verify that the test engineer is aware of all highly hazardous, reactive, or toxic components of the test material. The *test engineer* shall direct the test operator in proper safety procedures concerning these test materials.



- <u>4.4.2.</u> **Turn on** warning lights to indicate that hazardous testing is in progress. This places the test area in a *limited-access* condition.
- 4.4.3. **Ensure** that adequate supplies of LN<sub>2</sub> and LOX are available for testing purposes. *When needed*, **change out** the empty LOX dewars, and **cap** both the tap







and vent. Use hearing protection when moving or releasing LOX from dewars. Ensure that the flex line and diffuser tube do not hit the floor. Report as unusable any parts whose cleanliness is questionable; send these to the responsible cleaning agency.

- <u>4.4.4.</u> **Don** PPE (safety shoes, safety glasses, latex gloves).
- 4.4.5. **Assemble** the pre-cooling moat, and **install** it on the test stand as follows:

**CAUTION:** Verify proper operation of the safety catch before working in the plummet fall/impact path.

- a. Place the LN, moat atop the baseplate so that the hex head bolt holes are aligned and the open gap in the moat is facing the control room (north).
- b. Place the anvil into the LN<sub>2</sub> moat so that the bolts align with the hex head bolt holes in the bottom of the moat.
- c. **Place** the sample holder top plate on top of the anvil.
- d. **Insert** two bolts through the bolt holes.
- e. Tighten the two 9/16-in. hex head hold-down bolts, torquing each to 10 ft-lb (120 in.-lb) maximum.
- f. Screw in the striker pin guide, aligning the stud (northeast corner) and the striker pin guide swivel (southeast corner).
- **Place** the striker pin guide over the aligning stud and the guide swivel.
- h. **Place** the striker pin guide swivel spring onto the guide swivel, on top of the striker pin guide, and **screw** on the striker pin guide swivel nut.
- **Insert** the LN, feed-line nipple through the opening in the LN, moat. **Place** the nut and washer onto the nipple, and **tighten** the nut snugly by hand.
- **Place** the metal splash guard, flange side down, over the LN, nipple where it comes through the LN, moat.
- 4.4.6. Adjust the magnet housing assembly to the height needed to achieve the desired energy level (section 7.2, Table 1) by performing one of the following procedures:
  - 4.4.6.1. Procedure 1
    - 4.4.6.1.1. Load a blank and a striker pin in the test position.
    - 4.4.6.1.2. **Place** one end of the rod labeled with the appropriate heights on top of the striker pin.
    - 4.4.6.1.3. Ensure that the other end of the rod is in contact with the nose (lower end) of the plummet.
    - 4.4.6.1.4. **Tighten** the magnet housing in place.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 9 of 34

<u>4.4.6.1.5.</u> **Remove** the rod.

#### 4.4.6.2. Procedure 2

- 4.4.6.2.1. Load a blank and a striker pin in the test position.
- 4.4.6.2.2. In the control room, **use** the magnet switch on the control console to **release** the plummet so that it falls onto the safety catch.
- <u>4.4.6.2.3.</u> In the test cell, **support** the plummet, and **use** the foot switch to disengage the safety catch.
- <u>4.4.6.2.4.</u> Carefully lower the plummet until it is resting on top of the striker pin.
- 4.4.6.2.5. **Place** one end of the rod labeled with the proper heights on top of the plummet.
- <u>4.4.6.2.6.</u> **Ensure** that the other end of the rod is in contact with the bottom of the magnet housing.
- <u>4.4.6.2.7.</u> **Tighten** the magnet housing in place.
- <u>4.4.6.2.8.</u> **Remove** the rod.

**Note:** The proximity switch for the timer moves with the magnet housing assembly to the proper position for the test.



- 4.4.7. **Turn on** power for the control panel. **Ensure** the timer unit is plugged in. **Turn on** the power supply for the magnet.
- <u>4.4.8.</u> **Start** LN<sub>2</sub> flow to cool the base using either facility LN<sub>2</sub> or LN<sub>2</sub> dewars, as described below:

#### <u>4.4.8.1.</u> When using facility $LN_2$ :

- Ensure that the line connecting the LN<sub>2</sub> dewar is disconnected and capped.
- Open the LN<sub>2</sub> valve behind the test cell in Room 122. Don latex gloves, laboratory coat and safety glasses, and open the moat valve on the base to allow LN<sub>2</sub> to fill the moat and cool the base. As an option in warm weather, open the main LN<sub>2</sub> valve behind the test cell in Room 118 to speed the process of chilling the base. This step can be performed earlier in the procedure.



Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 10 of 34

Run LN<sub>2</sub> until the base is capable of holding LOX; allow approximately
45 minutes of chilling time, depending on the outdoor temperature. Pause
during testing to rechill the base as needed.

## 4.4.8.2. When running on LN, dewars:

- **Ensure** that the LN<sub>2</sub> valve behind the test cell in Room 122 is closed before uncapping the dewar line.
- Uncap the dewar line, and attach it to the liquid port of the LN, dewar.
- **Open** the LN<sub>2</sub> valve on the dewar. **Don** PPE; then **open** the moat valve on the base to allow LN<sub>2</sub> to fill the moat and cool the base.
- Run LN<sub>2</sub> until the base is capable of holding LOX; allow approximately
  45 minutes of chilling time, depending on the outdoor temperature. Pause
  during testing to rechill the base as needed.
- 4.4.9. While the base is cooling, **prepare** the striker pins and the sample freezing box. **Use** either of the two methods described below:

#### WARNING: Do not use glass dewars for transfer of cryogenic liquids.

### 4.4.9.1. Standard Chilling Method

- **Fill** an additional clean dewar with LOX for ladling into the sample tray and specimen cups.
- Place pins in a small clean dewar, and fill the dewar with enough LN<sub>2</sub> to cover pins with approximately 2 in. of liquid.
- **Ensure** that the sample freezing box is level. **Fill** the sample freezing box with LN<sub>2</sub>. **Place** the sample tray containing 20 prepared samples and a minimum of 5 blanks in the sample freezing box. **Wait** approximately 30 minutes, and **fill** the sample tray slowly with LOX until the sample cups are filled and the tray is filled to within 1/8 in. of its top.
- **Reject** any samples showing voids or holes. **Check** that the sample material has not separated from the bottom of the cup. Using tongs, **submerge** samples that separate. **Do not reject** cracked samples.
- **Ensure** that storage at this stage of the procedure does not exceed 4 hours and that sample cleanliness is maintained.

# 4.4.9.2. Alternate Chilling Method

- **Inspect** the striker pins for cleanliness and the striker-pin surface for smoothness. **Refer** to section 9.3, Required Tester Maintenance, *if there are any concerns*, and **inform** the test engineer. **Place** 25 clean, smooth pins in a small, clean, hand dewar.
- Partially fill the clean hand dewar containing striker pins with  $LN_2$ ; cover all pins with 2 in. of  $LN_2$ .





Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 11 of 34

- **Fill** an additional clean hand dewar with LOX for ladling into the specimen cups and the sample well.
- **Ensure** that the sample freezing box is level. **Remove** 20 prepared samples and a minimum of 5 blanks from the pan. **Place** these in the sample freezing box over the holes in the retainer plate. **Begin filling** the freezing box with LOX. **Fill** the freezing box and the cups to 1/8 in. from the top of the box.
- **Reject** any samples showing voids or holes. **Ensure** that sample material has not separated from the bottom of the cup. Cracked samples do not necessarily have to be rejected. **Inform** the test engineer of any problems.
- **Ensure** that storage at this stage of the procedure does not exceed 4 hours and that sample cleanliness is maintained.
- <u>4.4.10.</u> **Complete** the Pre-Test Checklist (section 7.2, Figure 1).
- 4.4.11. **Turn off** the moat valve. **Don** face shield, laboratory coat, and latex gloves; then **blow** ice and debris from the sample well with filtered GN<sub>2</sub>. **Blow out** all LN, to prevent dilution of LOX with LN<sub>2</sub>.
- <u>4.4.12.</u> **Check** the timing sensor before each drop to ensure proper positioning on the magnet base in relation to the plummet. The light on the sensor **shall be off** before the test begins.
- <u>4.4.13.</u> **Don** latex gloves. **Use** steel forceps to place a chilled specimen cup containing a chilled sample or blank into the base cavity. The first two cups **shall** be blanks to ensure the cleanliness of the system.
- <u>4.4.14.</u> Using steel forceps, **place** the striker pin on top of the blank or sample.

**Note:** *If aligning striker pin with gloved hands,* **do not touch** the impact surface of the pin.

- <u>4.4.15.</u> **Check** the striker-pin guide assembly to ensure that both guide pins are perpendicular to the guide plate, that the guide plate rotates freely with the spring, and that the guide plate is not bowed.
- <u>4.4.16.</u> **Slide** the striker-pin guide over the striker pin, and **ensure** that the striker pin is perpendicular to the sample.
- 4.4.17. Ladle LOX into the specimen cup, if the specimen is not covered with LOX.







Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 12 of 34

#### 4.5 Detailed Test Procedure

The *test operator* **shall perform** the following steps:

- 4.5.1. After completing system setup and sample loading, **ensure** there is no illumination other than the overhead light in the cell. **Proceed** to the control room, ensuring that the door between the test cell and control room is closed securely. After this step, *no access* to the test area shall be allowed.
- 4.5.2. At the control console, **reset** the **RED** velocity switch and the timer. **Ensure** that the timer reads **0.000 s** and is not running.
- 4.5.3. **Place** hand on the catch and magnet switches, and **ensure** an unobstructed view of the sample cup/striker pin configuration.
- 4.5.4. Turn off all lights in the test cell and control room area. (The light switches are to the right of the observation window.) Release the safety latch and magnet, in that order, by pushing up on the switches. Observe and listen to the drop. Count the number of rebounds, and note if a flash or audible report occurs on impact, first rebound, second rebound, etc.
- 4.5.5. **Turn on** all lights. On the test data sheet (section 7.2, Figure 2), **record** the drop time, any observations (visible flash, audible report, etc.), and when these observations occurred (first impact, first rebound, second rebound, etc.). The drop time **shall be** within 3% error, in accordance with Table 1, section 7.2.

**CAUTION:** Verify proper operation of the safety latch before working in the plummet fall/impact path.

- 4.5.6. Limited access to the test area is allowed at this point. **Don** latex gloves, clean laboratory coat, safety glasses, and any other personal protective equipment identified by the sample MSDS. **Re-enter** the test cell, **reposition** the plummet, and **verify** engagement of the safety latch. *If there is a noticeable odor*; **note** this on the test report.
- 4.5.7. **Lift** the striker-pin guide, **remove** the striker pin, and **observe** the striking surface. Then, **use** forceps to remove the specimen cup. **Save** the sample as tested, so that after the sample has warmed, the sample, specimen cup, inserts, and any substrates **shall** be examined for signs of obvious charring.
- 4.5.8. **Repeat** steps 4.4.11 through 4.4.17 and steps 4.5.1 through 4.5.7 for samples 2 through 20. **Test** a blank after every 5 samples to ensure that the test system is still clean. *If any samples react*, **test two blanks** after the reacted sample before continuing with the remaining samples. Any time a reaction occurs when dropping on a blank, **cease** all testing, and **blow** debris from the sample well with filtered GN<sub>2</sub>. **Resume** testing with two blank tests. *If a reaction occurs again*





Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 13 of 34

on a blank, cease testing, thaw the equipment, clean the base and plummet, **reassemble** the tester, and **resume** testing with two blanks.

- 4.5.9. On the test data sheet, **enter** observations of charred or burned defrosted test items.
- 4.5.10. **Do not reuse** inserts and insert cups. If the sample reacted, **photograph** the sample, save the insert, and package it with the sample tested. If the sample did not react, discard the insert. Save and package all insert cups.
- 4.5.11. If another test with a different material is to be performed, change out the base.
- 4.5.12. The test plan shall define the energy thresholding procedure. The following describes a general energy thresholding procedure:

Start testing at an impact energy of 72 ft-lb unless otherwise specified. Test in decrements of impact energies (drop heights) as listed in section 7.2, Table 1. If no reactions occur in 20 drops or if 1 reaction occurs in 60 drops, testing is complete. Upon 2 reactions, **drop** to the next energy level as noted in section 7.2, Table 1. Continue to drop the energy level until no reactions are observed in 20 drops or 1 reaction is observed in 60 drops; When this condition is met, testing is completed.

If the sample is depleted before thresholding is complete, **note** this on the test plan, and **discontinue** testing. The test engineer, in consultation with NASA, shall decide all questionable reactions. **Report** any violent reactions or anomalies with testing to the test engineer, who, with NASA's approval, shall decide whether to continue testing.

#### 4.6 Shutdown Procedure

The *test operator* **shall perform** the following steps:

- 4.6.1. **Don** latex gloves, clean laboratory coat, safety shoes, goggles, and other personal protective equipment identified by the LN, and LOX MSDSs. Close the LN<sub>2</sub> valve at the source, and then **close** the LN<sub>2</sub> valve at the moat. Closing the valves in this order ensures that no LN<sub>2</sub> remains trapped in the line between the valves.







- 4.6.2. **Place** a fan so that it blows on the base. **Turn on** the dehumidifiers.
- 4.6.3. When the stand has thawed sufficiently, **lift** the metal plate, and **check** for LN<sub>2</sub>/LOX. If none is visible, **disassemble** the base by performing section 4.4.5, steps a through j, in reverse.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 14 of 34

- <u>4.6.4.</u> **Turn off** the facility warning beacon to indicate that *full access* to the test cell is allowed at this point.
- 4.6.5. Cover the LOX diffuser on the LOX dewar with aluminum foil.
- 4.6.6. As soon as practical, **perform** data recording and post-test photography in accordance with section 4.7.
- 4.6.7. **Secure** the test area, *i.e.*, **close** and **lock** all doors leading from the test area to the outside of the facility. **Turn off** all lights in the test area.

#### 4.7 Data Recording and Reduction; Post-Test Photography

The *test operator* **shall perform** the following steps:

- 4.7.1. Complete the test report data sheets (section 7.2, Figure 2). Place these sheets and the completed Pre-Test Checklist (section 7.2, Figure 1) in the test folder.
- 4.7.2. Photograph reacted samples, charred or melted test equipment, or any anomalies. Document these in writing on the test report data sheet. Do not photograph post-test samples that did not react. Take photographs as close to the samples as possible. More than one sample or reaction per photograph is acceptable, if the details of reactions are visible. Refer to the *Photography Operating Guide* for procedures for taking photographs. Place three copies of each post-test photograph in the test folder before returning the folder to the engineer. Photographs shall be retained indefinitely.

**Note:** *If there are several reactions and samples are hard to handle,* representative photos may be taken and labeled as such. The test engineer, in consultation with NASA, **shall decide** whether to make representative photos on a case-by-case basis.

**Package** samples and inserts from reacted samples in clear photography slide sleeves. **Label** the protector with the test request number, NASA-STD-6001 test type, temperature at which the test was performed, reactions per number of samples tested, pressure at which the test was performed, and the date. **Use** a red pen to label reacted samples; **use** a black pen to label all other samples. **Identify** the sample by sample number. **Return** the samples with the test folder to the test engineer for evaluation. The *test engineer* **shall return** samples to the *sample preparation technician* who **shall store** them for future reference.



Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 15 of 34

# 5.0 Notes

Custodians for EM10-OWI-CHM-032		
Master List and Document Control	EM10 Management Support Assistant	
Alternate Document Control	EM10 ISO Representative	
Records	Materials Test Branch ISO Representative	
Calibration	Materials Test Branch Calibration Contact	
Memoranda	Materials Test Branch ISO Representative	

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 16 of 34

# **6.0 Safety Precautions and Warning Notes**

#### Warning

Death, severe personal injury, or loss of major equipment may result if maintenance or operating procedures, techniques, restrictions, etc., are not followed exactly.

#### **6.1 Hazards**

**Safety shall have precedence over all activities.** Because of the nature of testing materials in a LOX environment, the testing system involves several hazards to the operator and facility. Personnel conducting ambient pressure LOX impact testing may be exposed to the following hazards:

- Handling and pouring of cryogenic fluids
- Exposure to an oxygen-deficient environment
- Flammability and health risks from volatile cleaning solvents
- Handling and moving of heavy parts of the test apparatus on a regular basis
- Potential of explosion and hazardous fume by-products from burning materials in a pure oxygen or oxygen-enriched environment
- Electrical shock related to operation of equipment
- Potential touch temperature risks from cryogenic LOX tests when handling test equipment
- Pressurized systems with nitrogen
- Potential of impact plummet disengaging during assembly or disassembly of the impact chamber, creating risks of dismemberment, dislocation, or fracture.

### **6.2 Safety Precautions**

- 6.2.1. Personnel shall **plan** test setup, testing, and shutdown so that at least one test operator is in the test area and one other person is in Building 4623 during normal business hours. After normal business hours and on weekends, a test engineer shall be in Building 4623 during all test activities. **No more than five personnel** shall be in the test area at any given time. Operation of tests shall comply with EM10-OWI-CHM-050, *Building 4623 Guidelines for Test Operations*.
- <u>6.2.2.</u> When personnel are working with the tester, the controls shall not be **operated**. A sign warning that personnel are working in the test cell shall be placed on the control console.
- <u>6.2.3.</u> Personnel shall **refer to the MSDS** for information on personal protective equipment for materials being handled (sample materials, solvents used, gaseous nitrogen, liquid nitrogen, and liquid oxygen) and **shall wear** safety apparel appropriate for test specimens and conditions:
  - <u>Safety shoes</u> when there is a danger of foot injuries from falling or rolling objects, objects piercing the sole of the shoe, or when feet may be exposed to an electrical hazard









Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 17 of 34

- <u>Clean laboratory coats</u> when working with enriched oxygen or other oxidizers, combustion by-products, compressed gases, or flammable solvents
- <u>Chemically resistant goggles and gloves</u> while cleaning test equipment and while working with solvents
- <u>Clean thermal gloves and goggles</u> when pouring, handling, or transferring cryogenic fluids. Hydrocarbon residue can contaminate equipment and affect oxygen compatibility test results or damage the tester.
- Face shield when blowing LN, out of the base
- A respirator when working with solvents in closed or poorly ventilated spaces.
   Note that the appropriate respirator shall be worn as indicated on the MSDS.
   Cartridge respirators are only good for the constituents listed on the filtration cartridge and for dust particle filtration. Personnel shall be qualified to use the respirator, and the respirator shall be supplied by MSFC.
- Hearing protection during testing
- <u>Safety glasses</u> at all times while in the test cells.
- <u>6.2.4.</u> Serious tissue damage can occur on exposure to cryogenic LOX or LN<sub>2</sub>, cold vapors, or cold equipment. *If injury occurs*, personnel **shall call 911** and **ask** for medical assistance. Bystanders can (**but are not required to**) do the following:



- If it is safe to do so, remove the person from the source of cold.
- In the event of full-body cryogenic exposure and if it is safe to do so, remove the person from the exposure atmosphere, and keep the person's airway open. Loosely wrap the person in a blanket until medical personnel arrive.
- **Do not remove** frozen gloves, shoes, or clothing.
- **Do not massage** affected part(s).
- **Do not expose** affected part(s) to temperatures above 112 °F (45 °C).
- **Do not apply** ice, snow, or ointments to affected part(s).
- <u>6.2.5.</u> **Smoking is not permitted** in Building 4623. The test area is generally an oxygen-enriched environment. Open flame or other high-temperature sources are not permissible in the testing area while enriched-oxygen conditions exist. Personnel **shall not smoke or expose clothing** to an open flame for 30 minutes after handling liquid or gaseous oxygen.



- <u>6.2.6.</u> The building warning system **shall be activated** for the duration of all testing, including pre- and post-test procedures. Personnel **shall evacuate** the test area immediately when the oxygen alarm sounds and lights flash.
- <u>6.2.7</u>. *In case of an uncontrollable LOX fire*, personnel **shall not try to extinguish the fire but shall evacuate** the area immediately and shall **call 911** to notify the fire department.

Marshall Space Flight Center Organizational Work Instruction				
EM10				
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A		
NASA-STD-6001: Test 13A				

- <u>6.2.8.</u> Nothing shall be stored in the test cell area other than parts or components of the testing apparatus that are designated as spare parts and the tools necessary for routine equipment maintenance. All other materials shall be removed from the test area. Any spare parts that will be exposed to an enriched-oxygen environment shall be placed in the secured inventory area.
- <u>6.2.9.</u> Personnel shall **ensure** that the ventilation system fan is on continuously during testing, including pre-test and post-test activities, to bring in fresh air and remove fumes and other combustion by-products from the test cell.
- <u>6.2.10.</u> All testing shall be performed **remotely**. No one shall be allowed in the test cell during a test. Only the test operator(s) shall be allowed in the test cell during pre-test and post-test activities. Sample technicians are allowed to prepare liquid samples for testing. Other people shall be allowed in the test cell before and after testing as allowed by the test engineer.



- <u>6.2.11</u>. Personnel shall **ensure** all electrical components, wiring, *etc.*, are in good condition and properly connected and grounded and shall **use caution** when operating any electrical equipment. Electrical devices shall not be operated when floors in the test cell are wet.
- <u>6.2.12.</u> Glass dewar flasks shall not be used for the transfer of cryogenic fluids. (**Refer** to step 4.4.9 of this work instruction.)
- 6.2.13. All equipment that will contact LOX shall be cleaned as described in section 9.3, Required Tester Maintenance, and section 4.3, Equipment Checkout, before the equipment is exposed to LOX.
- <u>6.2.14</u> When handling cylinders and dewars or making connections for compressed gases and/or liquids, personnel shall **refer** to *Working Safely with Compressed Gases and Cryogens* and *NSTC 313-Cryogenics Safety* and shall **comply** with the suggestions inside these presentations. (The test engineer has these resources.)
- <u>6.2.15.</u> The building warning lights **shall be checked** daily for proper operation.

#### 6.3 Special Hazards Associated with Compressed Gases and Liquids

- <u>6.3.1.</u> All operations involving compressed gases and liquids shall be conducted with at least 2 people, in visual contact, in the facility.
- <u>6.3.2.</u> All operating personnel shall be instructed on the nature of hazards associated with compressed gases and liquids.
- <u>6.3.3.</u> Before removal of any component of the system for servicing, the *operator* **shall secure and inspect** the system to ensure that no unsafe condition exists.

Marshall Space Flight Center Organizational Work Instruction					
EM10					
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A					
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 19 of 34			

<u>6.3.4.</u> Personnel shall perform continuous monitoring, *e.g.*, check operating pressures, look for leaks, listen for unusual noises, during all operations. Personnel shall ensure that oxygen leak levels are adequate throughout operations.

#### **6.4 Emergency Shutdown Procedure**

The tester does not have to be shut down to be considered safe in an emergency situation.

#### 6.5 Accident Reporting

- <u>6.5.1.</u> From a safe location, the *test operator* **shall call 911 immediately** and **notify** the EM10 Branch Chief.
- <u>6.5.2.</u> From a safe location, the EM10 Branch Chief **shall immediately report** the accident to the NASA Safety Monitor and the appropriate supervisor(s).

#### **6.6 Emergency Response Plan**

Emergency procedures and plans for Building 4623 are incorporated into this OWI and are stated in MPR 1040.3, *MSFC Emergency Plan*, current revision. Plans shall be modified if operations change in a significant manner.

#### 6.7 Mishap Reporting

All mishaps occurring in Building 4623 **shall be reported** to the test engineer, who shall report the mishap to the Building 4623 Safety Monitor. An initial verbal report **shall be made** within 8 hours, followed by a written report within 3 days. The EM10 Branch Chief **shall prepare** a managerial report within 7 days. Both reports **shall be reviewed** by the test operator's supervisor and by the NASA Safety Monitor. The detail and extent of the mishap report **shall** depend on the nature and extent of the damage. *If personnel injury or equipment damage does occur*, the mishap report **shall** be completed in accordance with MWI 8621.1, *Close Call and Mishap Reporting and Investigation Program, current revision*.

Marshall Space Flight Center Organizational Work Instruction				
EM10				
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A				
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 20 of 34		

# 7.0 Attachments, Data, Reports, and Forms

#### 7.1 Attachment: Striker Pin Boxing Guide

This section describes the proper procedure for boxing striker pins for LOX cleaning. When two or more boxes of pins (~600) accumulate, they **shall** be sent for LOX cleaning, reducing the amount of batch testing. This boxing guide has been developed to allow:

- shipping boxes to be stacked without the bottom boxes being crushed
- ease of handling
- moving by hand truck
- quick count of pins.

The test operator **shall perform** the following steps:

- 7.1.1. **Mark** the number of pins in each rack with a permanent marker on the side of each rack.
- 7.1.2. **Pack** the racks in the special-purpose fiberboard shipping box, FSN 8815-00-117-8249 (a white file folder box with hand holes on the sides).
- <u>7.1.3.</u> **Place** two rows of five racks across the width of the box. **Orient** the racks so that the pins are on their sides.
- <u>7.1.4.</u> **Place** two racks flat on top of the two lower rows. This makes a maximum of 12 racks per box.
- 7.1.5. **Do not tape or seal** the boxes, since the pin count must be verified before the box is accepted at LOX cleaning.

#### 7.2 Forms

Figures 1, 2, and 3 contain typical forms used in support of Test 13A. Table 1 defines the standard thresholding procedure for Test 13A. Table 2 provides test information for other possible test parameters.

Marshall Space Flight Center Organizational Work Instruction					
EM10					
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A			
NASA-STD-6001: Test 13A					

# Ambient Pressure LOX Impact Tester Pre-Test Checklist

Figure 1. Typical Test 13A Pre-Test Checklist.

	Request No.		Initi Yes	al No
1.	NASA-STD-6001 reviewed?			
2.	OWI reviewed?			
3.	All safety equipment available and w	vorking properly?		
4.	Test equipment checked out per OW	/۱?		
5.	Sample MSDS read?			
6.	Test plan read and changes noted a by engineer?	nd approved		
7.	Adequate supplies of LOX and LN2'	?		
8.	Pins checked for surface roughness	/contamination?		
9.	Facility warning beacon activated?			
Test	Operator	Date		
1/05		EM	10-F-CHI	M-030

Note: Representative Data Sheet. Refer to Forms Master List for current version.

Marshall Space Flight Center Organizational Work Instruction					
EM10					
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A			
NASA-STD-6001: Test 13A					

Figure 2. Typical Test 13A Test Data Sheet.

TEST 13A TEST DATA SHEET					
Sample preparation information					
N	5.	5			
		Project:	_		
Manufacturer:		Batch: Lot:	_		
			_		
Substrate Material:	Su	bstrate Thickness (in.):	_		
		ecimen Cup Material:	_		
Cure Time: Cure T	emperature (F):	Cure Pressure (psia):	-		
Ave. weight (g)	<ul><li>Length (in.) _</li><li>Thickness (in.)</li></ul>	) Width (in.) Inserts Used: Y \_ N			
Tes	st operator i	nformation			
Drop Height (in.) kg	-m: En	vironment:	-		
	Test res	sults			
			_		
Sample Drop Flash* Audible* Time (s)	Odor Char**	Reaction*** Comments			
1		1			
2		2			
3		3			
4		4			
5		5			
blank		blank			
6		6			
7		7			
8		8			
9		9			
10		10			
blank		blank			
11		11	_		
12		12			
13		13	_		
14   15		14	_		
blank		15 blank			
16		16			
17		17	_		
18		18	_		
19		19	_		
20		20	_		
*S1 = 1st Impact; R1 = 1st Rebound;					
**S = Sample; C = Cup; P = Pin; I = In	sert	Excess Sample?: Y N N			
***Complete last. 0 = No; 1 = Yes		Quantity of Excess:			
Test Summary (reactions/test):	at ft₋lh:	•	-		
1/05	at it-ib	EM10-F-CHM-01	3		

Note: Representative Data Sheet. Refer to Forms Master List for current version.

Marshall Space Flight Center Organizational Work Instruction				
EM10				
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A				
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 23 of 34		

#### Figure 3. Typical Calibration Statement.

Calibration Statement: Categories IV	and V Equipment
Calibration is required before use per (Calibration before use for each test series and by the Using Line Organization Calibration Contacts: EM10/James Perkins,	d periodic testing n)
User Name: Equipment Description:	
(attach multiple components sheets if n	
Manufacturer:  ECN: Serial No.:  Date of Calibration:  Type of Software and Version:	
Listing of Standards Associated with Calibration	:
Are standards National Institute of Standards an Technology (NIST) traceable?	d Y N
Did calibration meet equipment manufacturer's specifications?	Y N
Calibration was performed by:	
Remarks:	

Note: Representative Data Sheet. Refer to Forms Master List for current version.

EM10-F-CHM-018

1/05

Marshall Space Flight Center Organizational Work Instruction					
EM10					
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A			
NASA-STD-6001: Test 13A					

Table 1. Standard Free-Fall Times for Ambient LOX Tester.

Impact Energy <sup>1</sup>		Drop Height	Theoretical	Actual Free-fall	+3% Error	
kg-m	ft-lb	Note	(in.)	Free-fall Time (s)	Time (s)	(s)
9.9	72.0	2	43.3	0.475	0.442	0.455
7.6	55.0	3	33.0	0.415	0.379	0.390
5.5	40.0	4	24.0	0.354	0.319	0.329
3.4	25.0	5	15.0	0.280	0.248	0.255
1.4	10.0		6.0	0.177	0.144	0.148

#### Notes:

- 1. If no reactions are observed at a specific energy level, testing shall be conducted at the energy level directly above the passing energy level, i.e., if no reactions are observed at 10 ft-lb, increase energy to 25 ft-lb, if the material has not been tested at 25 ft-lb.
- 2. If 2 reactions occur in >10 drops, lower energy to 55.0 ft-lb. If 2 reactions occur in <10 but >5 drops, lower to 40.0 ft-lb. If 2 reactions occur in <5 drops, lower to 25.0 ft-lb.
- 3. If 2 reactions occur in >10 drops, lower energy to 40.0 ft-lb. If 2 reactions occur in <10 but >5 drops, lower to 25.0 ft-lb. If 2 reactions occur in <5 drops, lower to 10 ft-lb.
- 4. If 2 reactions occur in >10 drops, lower energy to 25.0 ft-lb. If 2 reactions occur in <10 drops, lower to 10.0 ft-lb.
- 5. If 2 reactions occur in ≤60 drops, lower energy to 10.0 ft-lb.

Table 2. Other Free-Fall Times for Ambient LOX Tester.

Impact Energy		Drop Height	Theoretical	Actual Free-fall	+3% Error
kg-m	ft-lb	(in.)	Free-fall Time (s)	Time (s)	(s)
8.9	65.0	39.0	0.451	0.416	0.428
8.3	60.0	36.0	0.433	0.399	0.410
6.9	50.0	30.0	0.395	0.360	0.371
6.2	45.0	27.0	0.375	0.340	0.350
4.8	35.0	21.0	0.331	0.296	0.305
4.2	30.0	18.0	0.306	0.274	0.282
2.8	20.0	12.0	0.250	0.218	0.224
2.1	15.0	9.0	0.217	0.184	0.189

Marshall Space Flight Center Organizational Work Instruction			
EM10			
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A			
NASA-STD-6001: Test 13A			

8.0 Records

Records for Ambient LOX Impact testing shall consist of memoranda that contain test results and that are stored electronically in the Materials and Processes Technical Information System (MAPTIS).

#### 8.1 Memoranda

Memoranda containing test results shall be retained indefinitely by EM10. These memoranda shall be stored electronically in the MAPTIS database and shall be accessible by test request number or memorandum number.

#### 8.2 Calibration Records

Not applicable to Test 13A.

#### 8.3 Maintenance of Records

- 8.3.1. Memoranda less than 10 years old shall be maintained in ready-access files in MAPTIS; memoranda 10 years old or older shall be automatically transferred to historical files.
- 8.3.2. Calibration records shall be maintained on site for a minimum of 10 years, filed and indexed by test request number. These shall be stored in a manner that will protect them, *e.g.*, in a test folder stored in a metal file cabinet. After 10 years, calibration records shall be transferred to historical files.
- <u>8.3.3.</u> The original test records shall be saved for a minimum of 5 years.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page <b>26</b> of <b>34</b>

# 9.0 Tools, Equipment, and Materials

#### 9.1 Standard Configuration of Tester

The ambient-pressure LOX mechanical impact tester (Figure 4) is fastened to a 2-ft x 2-ft x 2-ft concrete pedestal. The tester itself is approximately 5.5 ft high, with a 1-in. thick stainless-steel baseplate attached. A 5.5-in. x 6.6-in. moat, 3.5 in. deep, is bolted to the baseplate. A 20-lb ( $\pm 0.05$ -lb) drop-weight plummet is held between three upright columns and rides these columns on ball bearings. The bearings are attached to triangular webs at the top and bottom of the plummet body and bear upon the grooves in the upright columns. An electromagnet holds the plummet at the desired drop height. The magnet housing assembly can be moved vertically on a fourth column to a maximum height of 48 in. above the sample to obtain the required impact energy to the sample.

A pre-chilled specimen cup (Figure 5) containing one pre-chilled sample configured as described in EM10-OWI-CHM-042, *Test Sample Preparation for Testing in Building 4623*, is placed in the pre-chilled well, and LOX is ladled into the cup until the cup is filled. (Figures 6 and 7 show specimen hardware used in ambient impact testing. Refer to the Definitions section, page 3, for more information on these items.) The striker pin (Figure 8) is placed on top of the sample, and the striker pin guide is placed over the pin (Figure 9). The plummet is dropped onto the striker pin, which transfers energy to the sample in the specimen cup.

The hardware standard configuration is controlled by the *Ambient Mechanical Impact Tester Configuration Control Book*, which is controlled by NASA.

#### 9.2 Procedure for Deviations

Deviations to the baselined tester configuration **shall have** NASA written approval. It is the responsibility of the test engineer to obtain the written approval. After written approval is received, the change shall be added to the *Ambient Mechanical Impact Tester Configuration Control Book*.

#### 9.3 Required Tester Maintenance

The standard maintenance program for this test chamber and related control equipment is divided into weekly and as-required service. In addition, the program involves a maintenance log, calibration, and a required spare parts inventory. The test operator shall perform the following procedures:

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 27 of 34





Figure 4. Typical Ambient Pressure Mechanical Impact Tester.

Note: Illustration is representative. Actual appearance may vary.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		
NASA-STD-6001: Test 13A		

Figure 5. One-piece Specimen Cup.

Figure 6. Insert Cup.







Figure 7. Insert/Disk.

Figure 8. Striker Pin.







Figure 9. Typical Sample Test Configuration.



Note: Illustrations are representative. Actual appearances may vary.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		Revision No.: A
NASA-STD-6001: Test 13A		

<u>9.3.1.</u> Weekly Maintenance. **Check** the dewar fill tube weekly to ensure there are no obstructions or kinks in the fill tube that would hinder proper operation. **Ensure** the cleanliness of the test area.

#### 9.3.2. As-Required Maintenance

- <u>9.3.2.1.</u> Using a visual comparator and the drawing specifications, **check** the finish of the post-test striker pin and specimen cup impact surfaces. *If necessary*, **send** the striker pin to the Fabrication Group for resurfacing. Pins shall always remain in their original rack (tolerance range), unless they are damaged or out of specification and, therefore, shall be discarded. (MSFC specifications are that the minimum distance from the top of the shoulder to the striking surface shall be 1.315 in.) **Do not replace** a discarded pin with a pin not in the tolerance range (0.005 in.) of the rack.
- 9.3.2.2. Place post-test pins in their original racks, bag them, and tape the bags closed. (See section 7.1 for striker pin boxing guide.) Count and label the number of pins on the bag with a permanent marker. Prepare a minimum batch of 600 pins as needed to go to the Cleaning Facility. Pins shall be cleaned to MSFC-SPEC-164B, Specification for Cleanliness of Components for Use in Oxygen, Fuel, and Pneumatic Systems, Class 1, Level A.
- 9.3.2.3. Newly machined cups and inserts shall be cleaned to MSFC-SPEC-164B, Specification for Cleanliness of Components for Use in Oxygen, Fuel, and Pneumatic Systems, Class 1, Level A. **Prepare** a minimum batch of 500 newly machined cups to go to the Cleaning Facility. **Bag** cups in groups of 25.
- <u>9.3.3.</u> Maintenance Log. **Document** any maintenance to the test chamber or setup in the *Ambient-Pressure LOX Mechanical Impact Tester Maintenance Log* to provide a history of the tester. Any deviation to standard maintenance shall be documented by the test operator and approved on the maintenance log by the responsible test engineer.

#### 9.4 Calibration

None.

#### 9.5 Required Spare Parts Inventory

**Verify** that the spare parts listed in Table 3 are available at the beginning of each test request, so that the testing of a material can be completed as close to within 1 working day as possible.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A		

Table 3. Required Spare Parts Inventory for Test 13A.

<u>Part</u>	<u>Quantity</u>	<u>Drawing #/Description</u>
Bolt, Plummet Guide Bar		
Bottom Plate, Anvil	1	B-D-1 Detail #1
Top Plate, Sample Holder	1	B-D-1 Detail #2
Moat, Pre-Cooling		
Bolt, Anvil Assembly Hold Down		
Stud, Striker Pin Guide Aligning		
Stud, Striker Pin Guide Swivel		
Nut, Striker Pin Guide Swivel		
Spring, Striker Pin Guide Swivel	2	1.5" x .470"
Guide, Striker Pin	5	B-A-6
Screw, Plummet Bearing Retainer		
Bearing, Plummet		
Spring, Striker Pin Guide	1	N/A

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		
NASA-STD-6001: Test 13A		

# **10.0 Personnel Training**

The nature of testing that occurs in Building 4623 is complex and involves potential hazards; therefore, all test operators shall complete the requirements for Category 1 Credentials before conducting any test, and all tester maintenance personnel shall complete the requirements for Category 2 Credentials.

- Category 1 Credentials qualify personnel to perform basic test operations.
- Category 2 Credentials qualify personnel to maintain and modify testing apparatus.

#### Category 1 Credentials - Basic Operations

To obtain Category 1 Credentials, the test operator shall complete training in following areas:

- High-Pressure Systems Safety
- Oxygen Compatibility
- General Safe Laboratory Practices
- Safe Handling of Cryogenic Fluids (LN, and LOX)
- Hazardous Waste Disposal.

#### Category 1 Credentialling also requires:

- Successful completion of an annual physical examination conducted by the medical facility at Marshall Space Flight Center (or equivalent), including a hearing exam
- A demonstration of knowledge of the test and equipment by the completion of two successful test sets under the supervision of the test engineer.
- A demonstration of knowledge of the OWI. Candidate test operators shall thoroughly read the test OWI and sign a statment confirming that they have read and understand the OWI. Each shall be issued a personal copy of the OWI.
- Passing of a written test covering the OWI. The test shall be administered by the test engineer.

A copy of the written test, along with the signed statement and the training record, shall constitute verification of credentials. Training records shall be kept on file as proof of training. These records shall include training expiration dates and required refresher courses.

Category 1 Credentials shall expire after a period of 2 years. After that time, recredentialling shall be required.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032		Revision No.: A
NASA-STD-6001: Test 13A		Page <b>32</b> of <b>34</b>

#### **Category 2 - Tester Maintenance and Modifications**

Personnel seeking **Category 2 Credentials** shall become qualified and credentialled through training classes approved by the candidate's supervisor or through training classes completed during previous employment. Training in the following areas shall be required:

- Compressed Gases and Working with Compressed Gas Lines and Fittings
- Basic Electrical Wiring.

This training shall be achieved through classes approved by the candidate's supervisor or through classes completed during previous employment.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING EM10-OWI-CHM-032 Revision No.: A		
NASA-STD-6001: Test 13A	Date: November 17, 2005	Page 33 of 34

This page is deliberately blank.

Marshall Space Flight Center Organizational Work Instruction		
EM10		
AMBIENT LOX IMPACT TESTING	EM10-OWI-CHM-032	Revision No.: A
NASA-STD-6001: Test 13A		

#