

Specifications for electrolytic polishing equipment:

NASA/LaRC has a requirement for electrolytic transmission electron microscope sample preparation and seeks a solution that shall meet the minimum following specifications:

- Ability to polish 2.3 mm and 3.0 mm diameter samples
- Automatic polishing voltage determination
- Ability to store polishing settings in local database
- Ability to pre-thin large specimens
- Ability to electrolytically punch specimens out at appropriate diameter
- Simultaneous thinning from both sides of the specimen

***Equipment shall be shipped and delivered by 9/30/13***